OLYMPUS®

VE-31/VE-32

Stereo Microscopes



Affordable, Easy-to-Use, Portable Microscopes


Suitable for a wide range of laboratory, research, assembly and inspection in factories, as well as for educational purposes, the Models VE-31 and VE-32 redefine standards for stereo microscopes. Although exceptionally compact and lightweight, the VE-31 and VE-32 offer all the functions required for most applicationsincluding interpupillary distance adjustment, diopter adjustment, focusing tension adjustment and a slider-mounted aperture diaphragm.


Because the optical system is a product of the outstanding technology that made Olympus famous, you are assured a sharp, clear image every time.

- 2-stage slider-mounted aperture diaphragm adjusts brightness. contrast and depth-of-focus.
- ·Removable stage plate finished in black-and-white to allow selection of background according to specimen type.
- ·Allows direct mounting of LSGA-3 Epi-Illuminator on the microscope body.




·Accepts PM-6 Photomicrographic Camera for photomicrography of specimens.





SPECIFICATIONS

Model number	VE-31	VE-32
Total magnification	20X	40X
Eyepiece	G10X-T, F.N. 20	G20X-T, F.N. 10
Objective	2X, fixed, N.A. 0.09, working distance: 62mm	
Interpupillary distance	5076mmm	46—72mm
Diopter adjustment	±5 diopters	
Focusing adjustment range	48mm, with focusing tension adjustment	
Stage plate	Black-and-white, 80mmø	





Specifications are subject to change without notice.

Olympus business areas







Medical and health-care



Imaging and information area







Industrial applications area

FOCUS ON LIFE











OLYMPUS OPTICAL CO.,LTD. OLYMPUS OPTICAL CO. (EUROPA) GMBH. OLYMPUS AMERICA INC. 5 Comando Cerder Drive, MeMile, NY 11747-3157, U.S.A. ECorporate Denter Disse, Mehres, N. 1. 12 results, Visco.

OLYMPUS SINGAPORE PTE LTD.

Bit K 211 Hardonson Road #13-50, Handerson Instatted Park, Singapore 159552, Singapore OLYMPUS OPTICAL CO.(UK.) LTD. OLYMPUS AUSTRALIA PTY. LTD.
Unit 4, 112-118 Teisners Road North Plyde, NSW2113, Australia